AESCULAP® **SO**.line® KERRISON Knochenstanzen

Such **Q**uality!

AESCULAP® SQ.line® KERRISON

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
Referenzen

- Manual cutting endurance test. Three samples of two worst-case SQ.line® KERRISONS and previous Aesculap KERRISONS each were used to cut lamb spine 25 000 times. All SQ.line® KERRISONS were perfectly functioning after the test cycles without repair and thereby showed up to three times increased durability.
- Jamming test. Five samples of three different SQ.line® KERRISONS and previous Aesculap KERRISONS were tested on lamb spine with 50 repeated punches. The SQ.line® KERRISON showed no jamming for the duration of the test.
- Destructive test. The required destructive force was measured for five samples of two worst-case SQ.line® KERRISONS and previous Aesculap KERRISONS each, biting on a metal plate. The SQ.line® KERRISONS withstood significantly greater forces on the cutting edge than the corresponding previous Aesculap KERRISONS.
- Milk lubrication test. Three samples of the SQ.line® KERRISONS were lubricated with an emulsion of water and surgical milk instrument lubricant. Afterwards, the samples were able to successfully punch a test card 50 times without jamming between the main part and the slider.
- Performance test aged. Five samples of the most delicate SQ.line® KERRISON bone punch with three corresponding previous Aesculap KERRISONS were compared. After aging the SQ.line® KERRISONS by 500 reprocessing cycles, they successfully passed all defined test criteria for the wearing parts. The SQ.line® KERRISON outperformed the previous Aesculap KERRISONS.
- Gloss value measurement. The gloss values of two new and five aged SQ.line® KERRISONS were measured on three defined areas from two different angles each. The results were compared to the values of uncoated Aesculap KERRISONS. The processing had no significant influence on the reflection behavior. The Noir® plus coating reduced reflections on the slider by 70% to 86%, on the main part by 47% to 63% and on the handle by 19% to 76% depending on the direction of measurement.
- Handle comfort evaluation report. The SQ.line® KERRISON handle was compared in terms of comfort during use to its predecessor and one of its competitor's handle. It is based on the visualization of prevailing surface pressure distributions during the actuation by means of pressure measuring foil. The SQ.line® KERRISON showed smaller areas of high surface pressure and a more homogeneous pressure distribution.
- Usability test with neurosurgeons. The usability test with neurosurgeons was performed to validate multiple customer requirements in terms of the performance and the ergonomics of SQ.line® KERRISON. The test was performed with 12 neurosurgeons within a simulated clinical environment. The performance and the ergonomic properties were rated positive.
- Usability test with CSSD and OR personnel. The usability test with CSSD and OR personnel was performed to validate multiple customer requirements in terms of the reprocessing and handling of the upwards cutting SQ.line® KERRISON. Ten people with an average CSSD experience of more than five years were interviewed. The suitability for reprocessing and the handling were rated positive.
- The SQ.line® KERRISON shake test was performed to investigate the resistance of the Noir® plus surface to scratches in comparison with the existing bone punches. For this purpose a sieve basket and a SQ.line® KERRISON rack, both containing bone punches, were fixed on a vibrating plate. It could be demonstrated, that the scratching of surface is strongly reduced with the use of the SQ.line® KERRISON racks.
- A usability test for the SQ.line® KERRISON racks was performed to validate the realization of multiple user needs and other requirements. Ten Aesculap employees who have knowledge in storage and reprocessing of surgical instruments and are familiar with general OR techniques and OR environment were interviewed. All users rated the usability of the racks positive.

B. Braun Austria GmbH | Otto Braun-Straße 3-5 | 2344 Maria Enzersdorf
Tel. +43 2236 46541-0 | bbraun.austria@bbraun.com | www.bbraun.at
 www.facebook.com/bbraunat  www.bbraun.info/youtube
 www.instagram.com/bbraun_at  www.bbraun.info/in

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany